

Product Information

Product SKU:	CAB17576	Gene ID:	8800	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Mouse

Additional Information

Observed MW:	28kDa	Conjugate:	Unconjugated
Calculated MW:	28kDa	Isotype:	IgG

Immunogen Information

Background: This gene is a member of the PEX11 family, which is composed of membrane elongation factors involved in regulation of peroxisome maintenance and proliferation. This gene product interacts with peroxisomal membrane protein 19 and may respond to outside stimuli to increase peroxisome abundance. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene.

Recommended Dilution:	WB,1:500 - 1:2000
Synonyms:	PMP28; hsPEX11p; PEX11alpha; PEX11-ALPHA; PEX11A
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 110-180 of human PEX11A (NP_003838.1).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.