

# PEX11A Rabbit Polyclonal Antibody



CAB17576

## Product Information

### Size:

20uL, 50uL, 100uL, 200uL

### Observed MW:

28kDa

### Calculated MW:

### Applications:

WB

### Reactivity:

Mouse

## Antibody Information

### Recommended dilutions:

WB 1:500 - 1:2000

### Source:

Rabbit

### Isotype:

IgG

### Purification:

Affinity purification

## Protein Background

This gene is a member of the PEX11 family, which is composed of membrane elongation factors involved in regulation of peroxisome maintenance and proliferation. This gene product interacts with peroxisomal membrane protein 19 and may respond to outside stimuli to increase peroxisome abundance. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene. [provided by RefSeq, Oct 2012]

## Immunogen information

### Gene ID:

8800

### Uniprot

O75192

### Synonyms:

PEX11-ALPHA; PMP28; hsPEX11p; PEX11A

### Immunogen:


Recombinant fusion protein containing a sequence corresponding to amino acids 110-180 of human PEX11A (NP\_003838.1).

### Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Product Images

---



Western blot - PEX11A Polyclonal Antibody (CAB17576)