
Product Information

Product SKU:	CAB17580	Gene ID:	8861	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	46kDa	Conjugate:	Unconjugated
Calculated MW:	47kDa	Isotype:	IgG

Immunogen Information

Background: Enables LIM domain binding activity; RNA polymerase II-specific DNA-binding transcription factor binding activity; and enzyme binding activity. Involved in negative regulation of transcription, DNA-templated and positive regulation of transcription by RNA polymerase II. Located in chromatin. Part of beta-catenin-TCF complex.

Recommended Dilution: WB,1:500 - 1:2000 IHC-P,1:50 - 1:100

Synonyms: NLI; CLIM2; LDB-1; CLIM-2; LDB1

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-50 of human LDB1 (NP_003884.1).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.