

**LIMD1 Rabbit Polyclonal Antibody**  
**CAB17585**



**Product Information**

<b>Product SKU:</b>	CAB17585	<b>Gene ID:</b>	8994	<b>Size:</b>	20uL, 100uL
<b>Clone No:</b>	-	<b>Host Species:</b>	Rabbit	<b>Reactivity:</b>	Human,Mouse

**Additional Information**

<b>Observed MW:</b>	72kDa	<b>Conjugate:</b>	Unconjugated
<b>Calculated MW:</b>	72kDa	<b>Isotype:</b>	IgG

**Immunogen Information**

<b>Background:</b>	Predicted to enable transcription corepressor activity. Involved in several processes, including negative regulation of hippo signaling; regulation of gene expression; and response to hypoxia. Acts upstream of or within P-body assembly and gene silencing by miRNA. Located in several cellular components, including P-body; adherens junction; and focal adhesion. Part of RISC complex.
<b>Recommended Dilution:</b>	WB,1:500 - 1:2000
<b>Synonyms:</b>	LIMD1
<b>Purification Method:</b>	Affinity purification
<b>Immunogen:</b>	Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of human LIMD1 (NP_055055.1).
<b>Storage:</b>	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.