

ITGBL1 Rabbit Polyclonal Antibody



CAB17592

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

54kDa

Calculated MW:

38kDa/43kDa/53kDa

Applications:

WB

Reactivity:

Human, Mouse, Rat

Protein Background

This gene encodes a beta integrin-related protein that is a member of the EGF-like protein family. The encoded protein contains integrin-like cysteine-rich repeats. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Nov 2012]

Immunogen information

Gene ID:

9358

Uniprot

O95965

Synonyms:

ITGBL1; OSCP; TIED

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Immunogen:

A synthetic peptide corresponding to a sequence within amino acids 200-300 of human ITGBL1 (NP_004782.1).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images

Western blot - ITGBL1 Polyclonal Antibody (CAB17592)

