

ARHGEF17 Rabbit Polyclonal Antibody



CAB17599

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

Refer to figures

Calculated MW:

Applications:

WB

Reactivity:

Human, Mouse

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

9828

Uniprot

Q96PE2

Synonyms:

P164RHOGEF; RHOGEF17; TEM4; p164-RhoGEF; ARHGEF17

Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 1150-1450 of human ARHGEF17 (NP_055601.2).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images
