

Product Information

Product SKU:	CAB17600	Gene ID:	9904	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human

Additional Information

Observed MW:	107kDa	Conjugate:	Unconjugated
Calculated MW:	107kDa	Isotype:	IgG

Immunogen Information

Background: This gene encodes a nucleolar protein that contains six RNA-binding motifs. The encoded protein may be involved in regulating ribosome biogenesis. Multiple alternatively spliced variants, encoding the same protein, have been identified.

Recommended Dilution:	WB,1:500 - 1:2000
Synonyms:	Mrd1; RBM19
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 759-960 of human RBM19 (NP_057280.2).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.