

LPCAT3 Rabbit Polyclonal Antibody



CAB17604

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

56kDa

Calculated MW:

Applications:

WB

Reactivity:

Mouse, Rat

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

10162

Uniprot

Q6P1A2

Synonyms:

C3F; LPCAT; LPLAT 5; LPSAT; MBOAT5; OACT5; nesy; LPCAT3


Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 160-370 of human LPCAT3 (NP_005759.4).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Western blot - LPCAT3 Polyclonal Antibody (CAB17604)