

ANGPTL7 Rabbit Polyclonal Antibody



CAB17606

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

Refer to figures

Calculated MW:

40kDa

Applications:

WB

Reactivity:

Human

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

10218

Uniprot

O43827

Synonyms:

ANGPTL7; AngX; CDT6; dJ647M16.1

Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 27-346 of human ANGPTL7 (NP_066969.1).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images
