

# NMUR1 Rabbit Polyclonal Antibody



**CAB17610**

---

## Product Information

### Size:

20uL, 50uL, 100uL, 200uL

### Observed MW:

53kDa

### Calculated MW:

47kDa

### Applications:

WB

### Reactivity:

Human, Mouse

## Antibody Information

### Recommended dilutions:

WB 1:500 - 1:2000

### Source:

Rabbit

### Isotype:

IgG

### Purification:

Affinity purification

## Protein Background

### Immunogen information

#### Gene ID:

10316

#### Uniprot

Q9HB89

#### Synonyms:

(FM-3); FM-3; FM3; GPC-R; GPR66; NMU1R; NMUR1

#### Immunogen:


Recombinant fusion protein containing a sequence corresponding to amino acids 360-426 of human NMUR1 (NP\_006047.3).

#### Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Product Images

---



Western blot - NMUR1 Polyclonal Antibody (CAB17610)