

Product Information

Product SKU:	CAB17615	Gene ID:	10570	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	62kDa	Conjugate:	Unconjugated
Calculated MW:	62kDa	Isotype:	IgG

Immunogen Information

Background: Enables filamin binding activity. Predicted to be involved in nervous system development. Predicted to be located in cytosol.

Recommended Dilution:	WB,1:500 - 1:2000
Synonyms:	CRMP3; DRP-4; ULIP4; CRMP3/DPYSL4
Purification Method:	Affinity purification
Immunogen:	A synthetic peptide corresponding to a sequence within amino acids 500-572 of human CRMP3/DPYSL4 (NP_006417.2).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.