

Product Information

Product SKU:	CAB17620	Gene ID:	10668	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human

Additional Information

Observed MW:	38kDa	Conjugate:	Unconjugated
Calculated MW:	38kDa	Isotype:	IgG

Immunogen Information

Background: Predicted to enable metal ion binding activity. Predicted to be involved in negative regulation of cell growth. Located in endoplasmic reticulum and nucleoplasm.

Recommended Dilution:	WB,1:500 - 1:2000
Synonyms:	CGR19; RNF197; CGRRF1
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 40-170 of human CGRRF1 (NP_006559.1).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.