

B4GAT1 Rabbit Polyclonal Antibody



CAB17624

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

Refer to figures

Calculated MW:

Applications:

WB

Reactivity:

Human, Mouse, Rat

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

This gene encodes a member of the beta-1, 3-N-acetylglucosaminyltransferase family. This enzyme is a type II transmembrane protein. It is essential for the synthesis of poly-N-acetyllactosamine, a determinant for the blood group i antigen. [provided by RefSeq, Jul 2008]

Immunogen information

Gene ID:

11041

Uniprot

O43505

Synonyms:

B3GN-T1; B3GNT1; B3GNT6; BETA3GNT1; MDDGA13; iGAT; iGNT; B4GAT1

Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 35-180 of human B4GAT1 (NP_006867.1).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images
