

BAZ2A Rabbit Polyclonal Antibody



CAB17630

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

280-310kDa

Calculated MW:

Applications:

WB IHC IF

Reactivity:

Human, Mouse, Rat

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000 IHC 1:50
- 1:100 IF 1:50 - 1:200

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

11176

Uniprot

Q9UIF9

Synonyms:

TIP5; WALp3; BAZ2A / TIP5; BAZ2A


Immunogen:

Recombinant fusion protein of Human BAZ2A.

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Western blot - BAZ2A Polyclonal Antibody (CAB17630)