

CAB17645

Product Information

Product SKU:	CAB17645	Gene ID:	22936	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Rat

Additional Information

Observed MW:	75kDa	Conjugate:	Unconjugated
Calculated MW:	72kDa	Isotype:	IgG

Immunogen Information

**Background:** Involved in snRNA transcription by RNA polymerase II. Located in nucleoplasm. Part of transcription elongation factor complex.

Recommended Dilution:	WB,1:500 - 1:1000
Synonyms:	MRCCAT1; ELL2
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 200-292 of human ELL2 (NP_036213.2).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.