

PCNX1 Rabbit Polyclonal Antibody



CAB17647

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

290kDa

Calculated MW:

Applications:

WB

Reactivity:

Human

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

This gene encodes an evolutionarily conserved transmembrane protein similar to the pecanex protein in Drosophila. The fly protein is a component of the Notch signaling pathway, which functions in several developmental processes. [provided by RefSeq, Jul 2016]

Immunogen information

Gene ID:

22990

Uniprot

Q96RV3

Synonyms:

PCNX; PCNXL1; pecanex; PCNX1


Immunogen:

A synthetic peptide of Human PCNX1.

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Western blot - PCNX1 Polyclonal Antibody (CAB17647)