

Product Information

Product SKU:	CAB17647	Gene ID:	22990	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human

Additional Information

Observed MW:	290kDa	Conjugate:	Unconjugated
Calculated MW:	259kDa	Isotype:	IgG

Immunogen Information

**Background:** This gene encodes an evolutionarily conserved transmembrane protein similar to the pecanex protein in Drosophila. The fly protein is a component of the Notch signaling pathway, which functions in several developmental processes.

Recommended Dilution:	WB,1:500 - 1:2000
Synonyms:	PCNX; PCNXL1; pecanex; PCNX1
Purification Method:	Affinity purification
Immunogen:	A synthetic peptide corresponding to a sequence within amino acids 2250-2341 of human PCNX1 (NP_055797.2).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.