

Product Information

Product SKU:	CAB17648	Gene ID:	22994	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse

Additional Information

Observed MW:	145kDa	Conjugate:	Unconjugated
Calculated MW:	122kDa	Isotype:	IgG

Immunogen Information

Background:	Enables protein homodimerization activity. Involved in several processes, including intraciliary transport involved in cilium assembly; protein localization to centrosome; and regulation of centrosome duplication. Located in several cellular components, including ciliary transition zone; intercellular bridge; and microtubule organizing center. Colocalizes with centrosome.
Recommended Dilution:	WB,1:500 - 1:2000
Synonyms:	AZ1; ZA1; AZI1; CEP131
Purification Method:	Affinity purification
Immunogen:	A synthetic peptide corresponding to a sequence within amino acids 1000-1083 of human CEP131 (NP_001306157.1).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.