

CEP131 Rabbit Polyclonal Antibody



CAB17648

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

145kDa

Calculated MW:

117kDa/121kDa/122kDa

Applications:

WB

Reactivity:

Human, Mouse

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

22994

Uniprot

Q9UPN4

Synonyms:

AZ1; AZI1; ZA1; CEP131


Immunogen:

A synthetic peptide of Human CEP131.

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Western blot - CEP131 Polyclonal Antibody (CAB17648)