

Product Information

Product SKU:	CAB17649	Gene ID:	22998	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human

Additional Information

Observed MW:	140kDa	Conjugate:	Unconjugated
Calculated MW:	122kDa	Isotype:	IgG

Immunogen Information

Background: Enables myosin II head/neck binding activity. Involved in several processes, including cytoplasmic actin-based contraction involved in cell motility; positive regulation of stress fiber assembly; and regulation of focal adhesion assembly. Located in stress fiber. Colocalizes with myosin II complex.

Recommended Dilution:	WB,1:500 - 1:2000
Synonyms:	LMO7B; LIMCH1A; LIMCH1
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 150-450 of human LIMCH1 (NP_001106189.1).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.