

# CLCC1 Rabbit Polyclonal Antibody



CAB17656

---

## Product Information

### Size:

20uL, 50uL, 100uL, 200uL

### Observed MW:

62kDa

### Calculated MW:

### Applications:

WB IHC

### Reactivity:

Human, Mouse

## Antibody Information

### Recommended dilutions:

WB 1:500 - 1:2000 IHC 1:50  
- 1:200

### Source:

Rabbit

### Isotype:

IgG

### Purification:

Affinity purification

## Protein Background

### Immunogen information

#### Gene ID:

23155

#### Uniprot

Q96S66

#### Synonyms:

MCLC; CLCC1

#### Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 302-501 of human CLCC1 (NP\_055942.1).

#### Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Product Images

---

Western blot - CLCC1 Rabbit pAb (CAB17656)

