

**CLCC1 Rabbit Polyclonal Antibody**

**CAB17656**



**Product Information**

<b>Product SKU:</b>	CAB17656	<b>Gene ID:</b>	23155	<b>Size:</b>	20uL, 100uL
<b>Clone No:</b>	-	<b>Host Species:</b>	Rabbit	<b>Reactivity:</b>	Human,Mouse

**Additional Information**

<b>Observed MW:</b>	62kDa	<b>Conjugate:</b>	Unconjugated
<b>Calculated MW:</b>	62kDa	<b>Isotype:</b>	IgG

**Immunogen Information**

**Background:** Predicted to enable chloride channel activity. Predicted to be involved in chloride transport. Located in endoplasmic reticulum and mitochondria-associated endoplasmic reticulum membrane. Implicated in retinitis pigmentosa 32.

<b>Recommended Dilution:</b>	WB,1:500 - 1:2000 IHC-P,1:50 - 1:200
<b>Synonyms:</b>	MCLC; RP32; CLCC1
<b>Purification Method:</b>	Affinity purification
<b>Immunogen:</b>	Recombinant fusion protein containing a sequence corresponding to amino acids 302-501 of human CLCC1 (NP_055942.1).
<b>Storage:</b>	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.