

FAM120A Rabbit Polyclonal Antibody



CAB17658

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

122kDa

Calculated MW:

68kDa/71kDa/117kDa/121kDa/125kDa

Applications:

WB

Reactivity:

Human

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

23196

Uniprot

Q9NZB2

Synonyms:

C9orf10; HBVPTPAP; OSSA; FAM120A; DNAPTP1; DNAPTP5


Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 1-140 of human FAM120A (NP_055427.2).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Western blot - FAM120A Polyclonal Antibody (CAB17658)