

Product Information

Product SKU:	CAB17658	Gene ID:	23196	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human

Additional Information

Observed MW:	122kDa	Conjugate:	Unconjugated
Calculated MW:	122kDa	Isotype:	IgG

Immunogen Information

Background: Enables RNA binding activity. Located in cytosol.

Recommended Dilution:	WB,1:500 - 1:2000
Synonyms:	OSSA; C9orf10; HBVPTPAP; FAM120A
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 1-140 of human FAM120A (NP_055427.2).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.