

PLCL2 Rabbit Polyclonal Antibody



CAB17660

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

126kDa

Calculated MW:

Applications:

WB

Reactivity:

Mouse

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

23228

Uniprot

Q9UPRO

Synonyms:

PLCE2; PLCL2

Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 130-300 of human PLCL2 (NP_055999.2).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images

Western blot - PLCL2 Polyclonal Antibody (CAB17660)

