

# HACL1 Rabbit Polyclonal Antibody



CAB17674

---

## Product Information

### Size:

20uL, 50uL, 100uL, 200uL

### Observed MW:

Refer to figures

### Calculated MW:

54kDa/56kDa/61kDa/63kDa

### Applications:

WB

### Reactivity:

Mouse

## Antibody Information

### Recommended dilutions:

WB 1:500 - 1:2000

### Source:

Rabbit

### Isotype:

IgG

### Purification:

Affinity purification

## Protein Background

### Immunogen information

#### Gene ID:

26061

#### Uniprot

Q9UJ83

#### Synonyms:

2-HPCL; HPCL; HPCL2; PHYH2; HACL1

#### Immunogen:

A synthetic peptide of Human HACL1.

#### Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Product Images

---