

FAM127B Rabbit Polyclonal Antibody



CAB17676

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

18kDa

Calculated MW:

Applications:

WB IF

Reactivity:

Human, Mouse

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000 IF 1:50 - 1:100

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

26071

Uniprot

Q9BWD3

Synonyms:

CXX1b; FAM127B; MAR8A; SIRH6


Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 1-113 of human FAM127B (NP_001071640.1).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Western blot - FAM127B Polyclonal Antibody (CAB17676)