

Product Information

Product SKU:	CAB17685	Gene ID:	27065	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	21kDa	Conjugate:	Unconjugated
Calculated MW:	-	Isotype:	IgG

Immunogen Information

Background: Predicted to enable clathrin light chain binding activity. Involved in apoptotic process. Located in endoplasmic reticulum.

Recommended Dilution:	WB,1:500 - 1:2000
Synonyms:	P21; D4S234; NEEP21; D4S234E; NSG1
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 1-80 of human NSG1 (NP_055207.1).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.