

NSG1 Rabbit Polyclonal Antibody



CAB17685

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

21kDa

Calculated MW:

Applications:

WB

Reactivity:

Human, Mouse, Rat

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

27065

Uniprot

Synonyms:

D4S234; D4S234E; NEEP21; P21; NSG1

Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 1-80 of human NSG1 (NP_055207.1).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images

Western blot - NSG1 Rabbit pAb (CAB17685)

