

FAM155B Rabbit Polyclonal Antibody



CAB17686

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

53kDa

Calculated MW:

Applications:

WB

Reactivity:

Mouse, Rat

Protein Background

This gene encodes a product belonging to a family of proteins with unknown function. The presence of two transmembrane domains suggests that this protein is a multi-pass membrane protein. [provided by RefSeq, Sep 2011]

Immunogen information

Gene ID:

27112

Uniprot

O75949

Synonyms:

CXorf63; TED; TMEM28; bB57D9.1; FAM155B

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Immunogen:

A synthetic peptide of Human FAM155B.

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images

Western blot - FAM155B Polyclonal Antibody (CAB17686)

