

# RNF115 Rabbit Polyclonal Antibody



**CAB17687**

---

## Product Information

### Size:

20uL, 50uL, 100uL, 200uL

### Observed MW:

34kDa

### Calculated MW:

### Applications:

WB

### Reactivity:

Rat

## Antibody Information

### Recommended dilutions:

WB 1:500 - 1:2000

### Source:

Rabbit

### Isotype:

IgG

### Purification:

Affinity purification

## Protein Background

### Immunogen information

#### Gene ID:

27246

#### Uniprot

Q9Y4L5

#### Synonyms:

BCA2; ZNF364; RNF115

#### Immunogen:


Recombinant fusion protein containing a sequence corresponding to amino acids 30-140 of human RNF115 (NP\_055270.1).

#### Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Product Images

---



Western blot - RNF115 Polyclonal Antibody (CAB17687)