

Product Information

Product SKU:	CAB17687	Gene ID:	27246	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Rat

Additional Information

Observed MW:	34kDa	Conjugate:	Unconjugated
Calculated MW:	34kDa	Isotype:	IgG

Immunogen Information

Background: Enables ubiquitin-protein transferase activity. Involved in negative regulation of epidermal growth factor receptor signaling pathway; protein autoubiquitination; and ubiquitin-dependent protein catabolic process via the multivesicular body sorting pathway. Predicted to be located in cytosol.

Recommended Dilution:	WB,1:500 - 1:2000
Synonyms:	BCA2; ZFP364; ZNF364; RABRING7; RNF115
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 30-140 of human RNF115 (NP_055270.1).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.