

# STRN4 Rabbit Polyclonal Antibody



CAB17694

---

## Product Information

### Size:

20uL, 50uL, 100uL, 200uL

### Observed MW:

98kDa

### Calculated MW:

### Applications:

WB

### Reactivity:

Rat

## Antibody Information

### Recommended dilutions:

WB 1:500 - 1:2000

### Source:

Rabbit

### Isotype:

IgG

### Purification:

Affinity purification

## Protein Background

### Immunogen information

#### Gene ID:

29888

#### Uniprot

Q9NRL3

#### Synonyms:

PPP2R6C; ZIN; zinedin; STRN4

#### Immunogen:


Recombinant fusion protein containing a sequence corresponding to amino acids 230-320 of human STRN4 (NP\_037535.2).

#### Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Product Images

---



Western blot - STRN4 Polyclonal Antibody (CAB17694)