

# EPN1 Rabbit Polyclonal Antibody



CAB17695

---

## Product Information

### Size:

20uL, 50uL, 100uL, 200uL

### Observed MW:

94kDa

### Calculated MW:

### Applications:

WB IHC

### Reactivity:

Human, Mouse, Rat

## Protein Background

This gene encodes a member of the epsin protein family. The encoded protein binds clathrin and is involved in the endocytosis of clathrin-coated vesicles. Loss of function of this gene is associated with reduced tumor growth and progression in certain cancer types. [provided by RefSeq, Mar 2016]

## Immunogen information

### Gene ID:

29924

### Uniprot

Q9Y6I3

### Synonyms:

EPN1

## Antibody Information

### Recommended dilutions:

WB 1:500 - 1:2000 IHC 1:50  
- 1:200

### Source:

Rabbit

### Isotype:

IgG

### Purification:

Affinity purification

### Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 130-180 of human EPN1 (NP\_037465.2).

### Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Product Images

---

Western blot - EPN1 Rabbit pAb (CAB17695)

