

MBIP Rabbit Polyclonal Antibody



CAB17707

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

39kDa

Calculated MW:

32kDa/34kDa/39kDa

Applications:

WB

Reactivity:

Rat

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

51562

Uniprot

Q9NS73

Synonyms:

MBIP


Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 1-344 of human MBIP (NP_057670.2).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Western blot - MBIP Polyclonal Antibody (CAB17707)