

Product Information

Product SKU:	CAB17707	Gene ID:	51562	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Rat

Additional Information

Observed MW:	39kDa	Conjugate:	Unconjugated
Calculated MW:	39kDa	Isotype:	IgG

Immunogen Information

Background: Enables identical protein binding activity and protein kinase inhibitor activity. Involved in histone H3 acetylation; positive regulation of JNK cascade; and positive regulation of gene expression. Located in cytosol; nucleolus; and nucleoplasm.

Recommended Dilution:	WB,1:500 - 1:2000
Synonyms:	MBIP
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 1-344 of human MBIP (NP_057670.2).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.