

# SETD4 Rabbit Polyclonal Antibody



CAB17708

---

## Product Information

### Size:

20uL, 50uL, 100uL, 200uL

### Observed MW:

50kDa

### Calculated MW:

### Applications:

WB

### Reactivity:

Human, Mouse

## Antibody Information

### Recommended dilutions:

WB 1:500 - 1:2000

### Source:

Rabbit

### Isotype:

IgG

### Purification:

Affinity purification

## Protein Background

### Immunogen information

#### Gene ID:

54093

#### Uniprot

Q9NVD3

#### Synonyms:

C21orf18; C21orf27; SETD4

#### Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 210-300 of human SETD4 (NP\_001007260).

#### Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Product Images

---

Western blot - SETD4 Polyclonal Antibody (CAB17708)

