

Product Information

Product SKU:	CAB17710	Gene ID:	54517	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	80kDa	Conjugate:	Unconjugated
Calculated MW:	75kDa	Isotype:	IgG

Immunogen Information

Background: Enables enzyme binding activity and pseudouridine synthase activity. Involved in several processes, including pseudouridine synthesis; regulation of hematopoietic stem cell differentiation; and regulation of mesoderm development. Located in nucleus.

Recommended Dilution:	WB,1:500 - 1:1000
Synonyms:	IDDABS; PUS7
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 352-661 of human PUS7 (NP_061915.2).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.