

PAK1IP1 Rabbit Polyclonal Antibody
CAB17713



Product Information

Product SKU:	CAB17713	Gene ID:	55003	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse

Additional Information

Observed MW:	44kDa	Conjugate:	Unconjugated
Calculated MW:	44kDa	Isotype:	IgG

Immunogen Information

Background: Involved in regulation of signal transduction by p53 class mediator and ribosomal large subunit biogenesis. Located in nucleolus.

Recommended Dilution:	WB,1:500 - 1:2000
Synonyms:	PIP1; MAK11; WDR84; hPIP1; bA421M1.5; PAK1IP1
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 1-392 of human PAK1IP1 (NP_060376.2).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.