

PAK1IP1 Rabbit Polyclonal Antibody



CAB17713

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

44kDa

Calculated MW:

43kDa

Applications:

WB

Reactivity:

Human, Mouse

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

55003

Uniprot

Q9NWT1

Synonyms:

MAK11; PIP1; WDR84; bA421M1.5; hPIP1; PAK1IP1


Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 1-392 of human PAK1IP1 (NP_060376.2).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Western blot - PAK1IP1 Polyclonal Antibody (CAB17713)