

Product Information

Product SKU:	CAB17714	Gene ID:	55079	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Rat

Additional Information

Observed MW:	49kDa	Conjugate:	Unconjugated
Calculated MW:	49kDa	Isotype:	IgG

Immunogen Information

Background: Predicted to enable transcription cis-regulatory region binding activity. Predicted to be involved in positive regulation of transcription, DNA-templated. Predicted to act upstream of or within several processes, including negative regulation of transcription by RNA polymerase II; nervous system development; and regulation of neuron differentiation. Predicted to be active in nucleus.

Recommended Dilution:	WB,1:500 - 1:2000
Synonyms:	FEZ; TOF; FEZL; FKSG36; ZFP312; ZNF312; FEZF2
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 1-280 of human FEZF2 (NP_060478.3).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.