

PI4K2B Rabbit Polyclonal Antibody



CAB17719

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

55kDa

Calculated MW:

Applications:

WB

Reactivity:

Human

Protein Background

This gene encodes a member of the type II PI4 kinase protein family. The encoded protein is primarily cytosolic and contributes to overall PI4-kinase activity along with other protein family members. This protein is involved in early T cell activation. [provided by RefSeq, Dec 2016]

Immunogen information

Gene ID:

55300

Uniprot

Q8TCG2

Synonyms:

PI4KIIB; PIK42B; PI4K2B

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification


Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 1-120 of human PI4K2B (NP_060793.2).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Western blot - PI4K2B Polyclonal Antibody (CAB17719)