

CAB17720

Product Information

Product SKU:	CAB17720	Gene ID:	55363	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse

Additional Information

Observed MW:	55kDa	Conjugate:	Unconjugated
Calculated MW:	55kDa	Isotype:	IgG

Immunogen Information

Background: Predicted to be involved in regulation of osteoblast differentiation. Located in nucleoplasm.

Recommended Dilution:	WB,1:500 - 1:2000
Synonyms:	NDR; EDAG; CT155; EDAG-1; HEMGN
Purification Method:	Affinity purification
Immunogen:	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human HEMGN (NP_932095.1).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.