

HEMGN Rabbit Polyclonal Antibody



CAB17720

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

55kDa

Calculated MW:

55kDa

Applications:

WB

Reactivity:

Human, Mouse

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

55363

Uniprot

Q9BXL5

Synonyms:

CT155; EDAG; EDAG-1; NDR; HEMGN; hemogen


Immunogen:

A synthetic peptide of Human HEMGN.

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Western blot - HEMGN Polyclonal Antibody (CAB17720)