

LYAR Rabbit Polyclonal Antibody



CAB17724

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

44kDa

Calculated MW:

43kDa

Applications:

WB IHC IF

Reactivity:

Human, Mouse, Rat

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000 IHC 1:50
- 1:200 IF 1:50 - 1:200

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

55646

Uniprot

Q9NX58

Synonyms:

LYAR; ZC2HC2; ZLYAR


Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 10-160 of human LYAR (NP_060286.1).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Western blot - LYAR Rabbit pAb (CAB17724)