

CAB17724

Product Information

Product SKU:	CAB17724	Gene ID:	55646	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	44kDa	Conjugate:	Unconjugated
Calculated MW:	44kDa	Isotype:	IgG

Immunogen Information

**Background:** Enables several functions, including DNA-binding transcription factor binding activity; identical protein binding activity; and transcription regulator inhibitor activity. Involved in several processes, including erythrocyte development; negative regulation of innate immune response; and regulation of transcription, DNA-templated. Located in nucleolus and nucleoplasm.

**Recommended Dilution:** WB,1:500 - 1:2000 IHC-P,1:50 - 1:200 IF/ICC,1:50 - 1:200

**Synonyms:** ZLYAR; ZC2HC2; LYAR

**Purification Method:** Affinity purification

**Immunogen:** Recombinant fusion protein containing a sequence corresponding to amino acids 10-160 of human LYAR (NP\_060286.1).

**Storage:** Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.