

UTP6 Rabbit Polyclonal Antibody



CAB17730

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

70kDa

Calculated MW:

Applications:

WB

Reactivity:

Human, Mouse, Rat

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

55813

Uniprot

Q9NYH9

Synonyms:

C17orf40; HCA66; UTP6


Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 1-220 of human UTP6 (NP_060898.2).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Western blot - UTP6 Polyclonal Antibody (CAB17730)