

DNAJC12 Rabbit Polyclonal Antibody



CAB17734

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

23kDa

Calculated MW:

12kDa/23kDa

Applications:

WB

Reactivity:

Human, Rat

Protein Background

This gene encodes a member of a subclass of the HSP40/DnaJ protein family. Members of this family of proteins are associated with complex assembly, protein folding, and export. Two transcript variants encoding distinct isoforms have been identified for this gene. [provided by RefSeq, Jul 2008]

Immunogen information

Gene ID:

56521

Uniprot

Q9UKB3

Synonyms:

DNAJC12; HPANBH4; JDP1

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification


Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 1-198 of human DNAJC12 (NP_068572.1).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Western blot - DNAJC12 Polyclonal Antibody (CAB17734)