

**CAB17736**

## **Product Information**

<b>Product SKU:</b>	CAB17736	<b>Gene ID:</b>	56902	<b>Size:</b>	20uL, 100uL
<b>Clone No:</b>	-	<b>Host Species:</b>	Rabbit	<b>Reactivity:</b>	Human,Rat

## **Additional Information**

<b>Observed MW:</b>	28kDa	<b>Conjugate:</b>	Unconjugated
<b>Calculated MW:</b>	28kDa	<b>Isotype:</b>	IgG

## **Immunogen Information**

**Background:** Enables RNA binding activity. Located in nucleolus and nucleoplasm.

<b>Recommended Dilution:</b>	WB,1:500 - 1:2000
<b>Synonyms:</b>	RRP20; KHRBP1; PNO1
<b>Purification Method:</b>	Affinity purification
<b>Immunogen:</b>	Recombinant fusion protein containing a sequence corresponding to amino acids 1-252 of human PNO1 (NP_064528.1).
<b>Storage:</b>	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.