

HHATL Rabbit Polyclonal Antibody



CAB17739

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

57kDa

Calculated MW:

Applications:

WB

Reactivity:

Human, Rat

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

57467

Uniprot

Q9HCP6

Synonyms:

C3orf3; GUP1; MBOAT3; MSTP002; OACT3; HHATL


Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 160-245 of human HHATL (NP_065758.3).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Western blot - HHATL Polyclonal Antibody (CAB17739)